

MPEG1 Antibody

Catalog No: #39997

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	MPEG1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	WB
Species Reactivity	Hu Rt Ms
Immunogen Type	Recombinant Protein
Immunogen Description	Synthesized peptide derived from human MPEG1
Target Name	MPEG1
Concentration	1.0mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	For extended storage, aliquot contents and freeze at -20 Λ C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Application Details

ELISA: 1:20000-1:80000 WB: 1:200-1:1000 IHC: 1:100-1:500

Images



Western blot analysis on different lysates using MPEG1 Antibody at 1/500 dilution.

Note: This product is for in vitro research use only