

MOXD1 Antibody

Catalog No: #39998



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description


Product Name	MOXD1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	ELISA WB IHC
Species Reactivity	Hu
Immunogen Type	Recombinant Protein
Immunogen Description	Rabbit polyclonal MOXD1 (1) antibody was raised against a recombinant human MOXD1 protein 196-368aa (BC018756).
Target Name	MOXD1
Concentration	0.54 mg/ml
Formulation	PBS, pH 7.4 with 0.02% sodium azide.
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Application Details

ELISA: 1:20000-1:80000 WB: 1:200-1:1000 IHC: 1:100-1:500

Images





Immunohistochemical analysis of paraffin-embedded MCF7 xenograft cytoplasmic staining using anti-MOXD1 Ab dilution at 1:100

Note: This product is for in vitro research use only