

MOCS3 Antibody


Catalog No: #40000

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MOCS3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Species Reactivity	Hu Rt Ms
Immunogen Type	Recombinant Protein
Immunogen Description	Rabbit polyclonal MOCS3 (1) antibody was raised against a recombinant human MOCS3 protein 15-175aa (BC015939).
Target Name	MOCS3
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% sodium azide.
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Images



Note: This product is for in vitro research use only