

MKRN3 Antibody

Catalog No: #40008

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MKRN3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen affinity column
Applications	ELISA WB IHC
Species Reactivity	Hu Rt Ms
Immunogen Type	Recombinant Protein
Immunogen Description	Rabbit polyclonal MKRN3 (1) antibody was raised against a recombinant human MKRN3 protein 278-443aa (BC044639).
Target Name	MKRN3
Concentration	1.0mg/ml
Formulation	Reconstitute in 100uL of sterile H ₂ O with 50% glycerol/0.01% thimerosal.
Storage	Store at +4°C for short term. Store at -20°C for long term. Avoid freeze/thaw cycle.


Application Details


WB 1:500-1:2000;

IHC 1:100-1:200;

ELISA 1:20000-1:80000

Images





Immunohistochemical staining of formalin-fixed paraffin-embedded fetal colon showing cytoplasmic staining with anti-MKRN3 antibody

Note: This product is for in vitro research use only