

CDX1 Antibody

Catalog No: #40057

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CDX1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total CDX1 protein.
Immunogen Type	Protein
Immunogen Description	Full length fusion protein
Target Name	CDX1
Accession No.	Swiss-Prot:P47902Gene Accssion:BC096252
Uniprot	P47902
GeneID	1044;
Concentration	0.8mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry:1:20-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human esophagus cancer tissue using #40057 at dilution 1/20.

Background

This gene is a member of the caudal-related homeobox transcription factor gene family. The encoded DNA-binding protein regulates intestine-specific gene expression and enterocyte differentiation. It has been shown to induce expression of the intestinal alkaline phosphatase gene, and inhibit beta-catenin/T-cell factor transcriptional activity.

Note: This product is for in vitro research use only