

SHPRH Antibody

Catalog No: #40104

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SHPRH Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total SHPRH protein.
Immunogen Type	Protein
Immunogen Description	Fusion protein corresponding to residues near the C terminal of human SNF2 histone linker PHD RING helicase, E3 ubiquitin protein ligase
Target Name	SHPRH
Other Names	bA54515.2
Accession No.	Swiss-Prot:Q149N8Gene Accssion:BC117686
Uniprot	Q149N8
GeneID	257218;
Concentration	2.4mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry:1:30-1:150

Images



Immunohistochemical analysis of paraffin-embedded Human lung cancer tissue using #40104 at dilution 1/45.

Background

SHPRH is a ubiquitously expressed protein that contains motifs characteristics of several DNA repair proteins, transcription factors, and helicases. SHPRH is a functional homolog of *S. cerevisiae* RAD5 (Unk et al., 2006 [PubMed 17108083]).

Note: This product is for in vitro research use only