

SRP54 Antibody

Catalog No: #40106

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	SRP54 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	WB IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous levels of total SRP54 protein.
Immunogen Type	Protein
Immunogen Description	Fusion protein of human signal recognition particle 54kDa
Target Name	SRP54
Accession No.	Swiss-Prot:P61011Gene Accssion:BC000652
Uniprot	P61011
GeneID	6729;
SDS-PAGE MW	56KD
Concentration	2.3mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:200-1:1000

Immunohistochemistry:1:25-1:100

Images




Gel: 8%SDS-PAGE

Lysate: 40ug Human bladder carcinoma tissuePrimary antibody: 1/500 dilution

Secondary antibody dilution: 1/8000

Exposure time: 2 minutes



Immunohistochemical analysis of paraffin-embedded Human prostate cancer tissue using #40106 at dilution 1/40.

Background

Binds to the signal sequence of presecretory protein when they emerge from the ribosomes and transfers them to TRAM (translocating chain-associating membrane protein).

Note: This product is for in vitro research use only