

THYN1 Antibody

Catalog No: #40147

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	THYN1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total THYN1 protein.
Immunogen Type	Protein
Immunogen Description	Full length fusion protein
Target Name	THYN1
Other Names	MY105; THY28; MDS012; HSPC144; THY28KD
Accession No.	Swiss-Prot:Q9P016Gene Accssion:BC006978
Uniprot	Q9P016
GeneID	29087;
SDS-PAGE MW	26KD
Concentration	1.5mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:200-1:1000

Images



Gel: 12%SDS-PAGE
 Lysate: 40ug Human bladder carcinoma tissue
 Primary antibody: 1/400 dilution
 Secondary antibody dilution: 1/8000
 Exposure time: 10 seconds

Background

This gene encodes a protein that is highly conserved among vertebrates and plant species and may be involved in the induction of apoptosis. Alternatively spliced transcript variants encoding different isoforms have been described.?

Note: This product is for in vitro research use only