

TRIM11 Antibody

Catalog No: #40148

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TRIM11 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total TRIM11 protein.
Immunogen Type	Protein
Immunogen Description	Fusion protein of human tripartite motif containing 11
Target Name	TRIM11
Other Names	BIA1; RNF92
Accession No.	Swiss-Prot:Q96F44Gene Accssion:BC011629
Uniprot	Q96F44
GeneID	81559;
Concentration	1.6mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry:1:25-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human breast cancer tissue using #40148 at dilution 1/30.

Background

The protein encoded by this gene is a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. This protein localizes to the nucleus and the cytoplasm. Its function has not been identified.?

Note: This product is for in vitro research use only