

HTATIP2 Antibody

Catalog No: #40151

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	HTATIP2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total HTATIP2 protein.
Immunogen Type	Protein
Immunogen Description	Fusion protein corresponding to residues near the C terminal of human HIV-1 Tat interactive protein 2, 30kDa
Target Name	HTATIP2
Other Names	CC3; TIP30; SDR44U1
Accession No.	Swiss-Prot:Q9BUP3Gene Accssion:BC002439
Uniprot	Q9BUP3
GeneID	10553;
SDS-PAGE MW	27KD
Concentration	2.1mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500-1:2000

Images



Gel: 8%SDS-PAGE
 Lysate: 40ug A431 cell
 Primary antibody: 1/800 dilution
 Secondary antibody dilution: 1/8000
 Exposure time: 5 seconds

Background

TIP30 (30 kDa HIV-1 TAT-interacting protein), also known as HTATIP2 (HIV-1 TAT-interactive protein 2) or CC3, is a ubiquitously expressed protein that functions as an oxidoreductase. It acts as a tumor suppressor participating in metabolic suppression, inhibition of angiogenesis and induction of

apoptosis.?

Note: This product is for in vitro research use only