

TMED3 Antibody

Catalog No: #40155

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TMED3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total TMED3 protein.
Immunogen Type	Protein
Immunogen Description	Fusion protein corresponding to residues near the C terminal of human transmembrane emp24 protein transport domain containing 3
Target Name	TMED3
Other Names	p26; P24B; C15orf22
Accession No.	Swiss-Prot:Q9Y3Q3Gene Accssion:BC000027
Uniprot	Q9Y3Q3
GeneID	23423;
Concentration	1.6mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry:1:25-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human breast cancer tissue using #40155 at dilution 1/30.

Background

Potential role in vesicular protein trafficking, mainly in the early secretory pathway. Contributes to the coupled localization of TMED2 and TMED10 in the cis-Golgi network.

Note: This product is for in vitro research use only