

TXNDC17 Antibody

Catalog No: #40172

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TXNDC17 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total TXNDC17 protein.
Immunogen Type	Protein
Immunogen Description	Full length fusion protein
Target Name	TXNDC17
Other Names	TRP14; TXNL5
Accession No.	Swiss-Prot:Q9BRA2Gene Accssion:BC006405
Uniprot	Q9BRA2
GeneID	84817;
Concentration	1.8mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry:1:25-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human breast cancer tissue using #40172 at dilution 1/30.

Background

Disulfide reductase. May participate in various redox reactions through the reversible oxidation of its active center dithiol to a disulfide and catalyze dithiol-disulfide exchange reactions. Modulates TNF-alpha signaling and NF-kappa-B activation. Has peroxidase activity and may contribute to the elimination of cellular hydrogen peroxide.

Note: This product is for in vitro research use only