

PLXNB2 Antibody

Catalog No: #40207

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PLXNB2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total PLXNB2 protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide corresponding to residues near the N terminal of human plexin B2
Target Name	PLXNB2
Other Names	MM1; PLEXB2; Nbla00445; dJ402G11.3
Accession No.	Swiss-Prot:O15031Gene Accssion:NP_036533
Uniprot	O15031
GeneID	23654;
Concentration	0.8mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry:1:25-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human thyroid cancer tissue using #40207 at dilution 1/20.



Immunohistochemical analysis of paraffin-embedded Human liver cancer tissue using #40207 at dilution 1/20.

Background

Members of the B class of plexins, such as PLXNB2 are transmembrane receptors that participate in axon guidance and cell migration in response to semaphorins (Perrot et al. (2002) [PubMed 12183458]).

Note: This product is for in vitro research use only