

EMC1 Antibody

Catalog No: #40208

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	EMC1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total EMC1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide corresponding to residues near the C terminal of human ER membrane protein complex subunit 1
Target Name	EMC1
Other Names	KIAA0090
Accession No.	Swiss-Prot:Q8N766Gene Accssion:NP_001258356
Uniprot	Q8N766
GeneID	23065;
SDS-PAGE MW	112KD
Concentration	1.1mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:200-1:1000

Images



Gel: 6%SDS-PAGE
 Lysate: 40ug Hepg2 cellPrimary antibody: 1/200 dilution
 Secondary antibody dilution: 1/8000
 Exposure time: 15 seconds

Background

This gene encodes a single-pass type I transmembrane protein, which is a subunit of the endoplasmic reticulum membrane protein complex (EMC).

Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene.

Note: This product is for in vitro research use only