

STAB1 Antibody

Catalog No: #40225

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	STAB1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total STAB1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide corresponding to residues near the C terminal of human stabilin 1
Target Name	STAB1
Other Names	FEX1; FEEL-1; FELE-1; STAB-1; CLEVER-1
Accession No.	Swiss-Prot:Q9NY15Gene Accssion:NP_055951
Uniprot	Q9NY15
GeneID	23166;
Concentration	1mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Colorectal tissue using #40225 at dilution 1/100



Immunohistochemical analysis of paraffin-embedded Human Thyroid tissue using #40225 at dilution 1/100

Background

This gene encodes a large, transmembrane receptor protein which may function in angiogenesis, lymphocyte homing, cell adhesion, or receptor scavenging. The protein contains 7 fasciclin, 16 epidermal growth factor (EGF)-like, and 2 laminin-type EGF-like domains as well as a C-type lectin-like hyaluronan-binding Link module. The protein is primarily expressed on sinusoidal endothelial cells of liver, spleen, and lymph node.?

Note: This product is for in vitro research use only