

TRIB1 Antibody

Catalog No: #40265

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TRIB1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total TRIB1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide corresponding to residues near the C terminal of human tribbles pseudokinase 1
Target Name	TRIB1
Other Names	C8FW; GIG2; TRB1; GIG-2; SKIP1; TRB-1
Accession No.	Swiss-Prot:Q96RU8 Gene Accssion:NP_079471
Uniprot	Q96RU8
GeneID	10221;
Concentration	1.4mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry:1:25-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human liver cancer tissue using #40265 at dilution 1/20.



Immunohistochemical analysis of paraffin-embedded Human breast cancer tissue using #40265 at dilution 1/20.

Background

Interacts with MAPK kinases and regulates activation of MAP kinases. May not display kinase activity.

Note: This product is for in vitro research use only