

USP9X Antibody

Catalog No: #40285

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	USP9X Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total USP9X protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide of human ubiquitin specific peptidase 9, X-linked
Target Name	USP9X
Other Names	FAF; FAM; DFFRX; MRX99
Accession No.	Swiss-Prot:Q93008 Gene Accssion:NP_001034679
Uniprot	Q93008
GeneID	8239;
Concentration	2.4mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry:1:25-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human esophagus cancer tissue using #40285 at dilution 1/35.



Immunohistochemical analysis of paraffin-embedded Human thyroid cancer tissue using #40285 at dilution 1/35.

Background

This gene is a member of the peptidase C19 family and encodes a protein that is similar to ubiquitin-specific proteases. Though this gene is located on the X chromosome, it escapes X-inactivation. Mutations in this gene have been associated with Turner syndrome. Alternate transcriptional splice variants, encoding different isoforms, have been characterized.

Note: This product is for in vitro research use only