

ANXA2R Antibody

Catalog No: #46929

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ANXA2R Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total ANXA2R protein.
Immunogen Type	peptide
Immunogen Description	Synthetic peptide of human ANXA2R
Target Name	ANXA2R
Other Names	AX2R; AXIIR; C5orf39
Accession No.	Swiss-Prot#:Q3ZCQ2NCBI Gene ID:389289Gene Accssion:NP_001014301
Uniprot	Q3ZCQ2
GeneID	389289;
Calculated MW	40 kDa
Concentration	1.9mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20C

Application Details

Western blotting:1:200-1000

Images



Gel: 12%SDS-PAGE

Lysate: 40 µg, Lane: Human small cell malignant tumors (lymphoma) tissue lysate

Primary antibody:ANXA2R Antibody at dilution 1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 10 minutes

Background

May act as a receptor for annexin II on marrow stromal cells to induce osteoclast formation.

Note: This product is for in vitro research use only