

AOAH Antibody

Catalog No: #46930

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	AOAH Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total AOAH protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human AOAH
Target Name	AOAH
Accession No.	Swiss-Prot#:P28039NCBI Gene ID:313Gene Accssion:BC025698
Uniprot	P28039
GeneID	313;
Concentration	0.8mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20C


Application Details

Immunofluorescence:1: 20-100

Images



The image is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using 46930(AOAH Antibody) at dilution 1/25. (Original magnification: ?00)



The image is immunohistochemistry of paraffin-embedded Human cervical cancer tissue using 46930(AOAH Antibody) at dilution 1/25. (Original magnification: 200)

Background

This locus encodes both the light and heavy subunits of acyloxyacyl hydrolase. The encoded enzyme catalyzes the hydrolysis of acyloxyacyl-linked fatty acyl chains from bacterial lipopolysaccharides, effectively detoxifying these molecules. The encoded protein may play a role in modulating host inflammatory response to gram-negative bacteria. Alternatively spliced transcript variants have been described.

Note: This product is for in vitro research use only