

CDK18 Antibody

Catalog No: #46946

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CDK18 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total CDK18 protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human CDK18
Target Name	CDK18
Other Names	PCTK3; PCTAIRE; PCTAIRE3
Accession No.	Swiss-Prot#:Q07002NCBI Gene ID:5129Gene Accssion:BC011526
Uniprot	Q07002
GeneID	5129;
Concentration	0.5mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20C

Application Details

Immunofluorescence:1: 20-100

Images



The image is immunohistochemistry of paraffin-embedded Human liver cancer tissue using 46946(CDK18 Antibody) at dilution 1/20. (Original magnification: 200)

Background

May play a role in signal transduction cascades in terminally differentiated cells.

Note: This product is for in vitro research use only