

CDX4 Antibody

Catalog No: #46947

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CDX4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total CDX4 protein.
Immunogen Type	peptide
Immunogen Description	Synthetic peptide of human CDX4
Target Name	CDX4
Accession No.	Swiss-Prot#:O14627 NCBI Gene ID:1046Gene Accssion:NP_005184
Uniprot	O14627
GeneID	1046;
Concentration	0.5mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20C


Application Details

Immunofluorescence:1: 20-100

Images



The image is immunohistochemistry of paraffin-embedded Human gastric cancer tissue using 46947(CDX4 Antibody) at dilution 1/20. (Original magnification: ?00)



The image is immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using 46947(CDX4 Antibody) at dilution 1/20. (Original magnification: ?00)

Background

This gene encodes a member of a small subfamily of homeobox containing transcription factors. The encoded protein may regulate homeobox gene expression during anteroposterior patterning and hematopoiesis.

Note: This product is for in vitro research use only