

CLUL1 Antibody

Catalog No: #46957

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CLUL1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total CLUL1 protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human CLUL1
Target Name	CLUL1
Other Names	RA337M
Accession No.	Swiss-Prot#:Q15846NCBI Gene ID:27098Gene Accssion:BC025381
Uniprot	Q15846
GeneID	27098;
Concentration	0.5mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20C

Application Details

Immunofluorescence:1: 20-100

Images



The image is immunohistochemistry of paraffin-embedded Human liver cancer tissue using 46957(CLUL1 Antibody) at dilution 1/20. (Original magnification: ?00)

Background

CLUL1 is a retina-specific clusterin-like protein which shows nearly 25% identity to clusterin at the protein level. It is expressed specifically in cone photoreceptor cells.

Note: This product is for in vitro research use only