

CNDP2 Antibody

Catalog No: #46958

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CNDP2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total CNDP2 protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human CNDP2
Target Name	CNDP2
Other Names	CN2; CPGL; PEPA; HsT2298; HEL-S-13
Accession No.	Swiss-Prot#:Q96KP4NCBI Gene ID:55748Gene Accssion:BC001375
Uniprot	Q96KP4
GeneID	55748;
Concentration	0.3mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20C


Application Details

Immunofluorescence:1: 25-100

Images



The image is immunohistochemistry of paraffin-embedded Human liver cancer tissue using 46958(CNDP2 Antibody) at dilution 1/20. (Original magnification: ?00)



The image is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using 46958(CNDP2 Antibody) at dilution 1/20. (Original magnification: ?00)

Background

CNDP2, also known as tissue carnosinase and peptidase A (EC 3.4.13.18), is a nonspecific dipeptidase rather than a selective carnosinase (Teufel et al., 2003 [PubMed 12473676]).

Note: This product is for in vitro research use only