

## CNOT6 Antibody

Catalog No: #46962

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	CNOT6 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total CNOT6 protein.
Immunogen Type	peptide
Immunogen Description	Synthetic peptide of human CNOT6
Target Name	CNOT6
Other Names	CCR4; Ccr4a
Accession No.	Swiss-Prot#:Q9ULM6 NCBI Gene ID:57472 Gene Accssion:NP_001290170
Uniprot	Q9ULM6
GeneID	:57472
Concentration	0.8mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20C


## Application Details

Immunofluorescence:1: 20-100

## Images



The image is immunohistochemistry of paraffin-embedded Human liver cancer tissue using 46962(CNOT6 Antibody) at dilution 1/20. (Original magnification: ?00)



The image is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using 46962(CNOT6 Antibody) at dilution 1/20. (Original magnification: ?00)

## Background

This gene encodes the catalytic component of the CCR4-NOT core transcriptional regulation complex. The encoded protein has a 3'-5' RNase activity and prefers polyadenylated substrates. The CCR4-NOT complex plays a role in many cellular processes, including miRNA-mediated repression, mRNA degradation, and transcriptional regulation.

Note: This product is for in vitro research use only