

COPS7A Antibody

Catalog No: #46971

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	COPS7A Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total COPS7A protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human COPS7A
Target Name	COPS7A
Other Names	CSN7; CSN7A; SGN7a
Accession No.	Swiss-Prot#:Q9UBW8NCBI Gene ID:50813Gene Accssion:BC011789
Uniprot	Q9UBW8
GeneID	50813;
Concentration	0.7mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20C


Application Details

Immunofluorescence:1: 25-100

Images



The image is immunohistochemistry of paraffin-embedded Human gastric cancer tissue using 46971(COPS7A Antibody) at dilution 1/30. (Original magnification: ?00)



The image is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using 46971(COPS7A Antibody) at dilution 1/30. (Original magnification: ?00)

Background

This gene encodes a component of the COP9 signalosome, an evolutionarily conserved multi-subunit protease that regulates the activity of the ubiquitin conjugation pathway. Alternatively spliced transcript variants that encode the same protein have been described.

Note: This product is for in vitro research use only