

CPSF2 Antibody

Catalog No: #46984

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

| | |
|-----------------------|---|
| Product Name | CPSF2 Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | Antigen affinity purification |
| Applications | IHC |
| Species Reactivity | Hu |
| Specificity | The antibody detects endogenous levels of total CPSF2 protein. |
| Immunogen Type | peptide |
| Immunogen Description | Synthetic peptide of human CPSF2 |
| Target Name | CPSF2 |
| Other Names | CPSF100 |
| Accession No. | Swiss-Prot#:Q9P2I0 NCBI Gene ID:53981 Gene Accssion:NP_059133 |
| Uniprot | Q9P2I0 |
| GeneID | 53981; |
| Concentration | 0.8mg/ml |
| Formulation | Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol. |
| Storage | Store at -20C |

Application Details

Immunofluorescence:1: 25-100

Images



The image is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using 46984(CPSF2 Antibody) at dilution 1/20. (Original magnification: ?00)

Background

Component of the cleavage and polyadenylation specificity factor (CPSF) complex that play a key role in pre-mRNA 3'-end formation, recognizing the AAUAAA signal sequence and interacting with poly(A) polymerase and other factors to bring about cleavage and poly(A) addition. Involved in the histone 3' end pre-mRNA processing.

Note: This product is for in vitro research use only