

CTR9? Antibody

Catalog No: #46994

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	CTR9? Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total CTR9? protein.
Immunogen Type	peptide
Immunogen Description	Synthetic peptide of human CTR9?
Target Name	CTR9?
Other Names	TSBP; p150; SH2BP1; p150TSP
Accession No.	Swiss-Prot#:Q6PD62NCBI Gene ID:9646Gene Accssion:NP_055448
Uniprot	Q6PD62
GeneID	9646;
Calculated MW	134 kDa
Concentration	1mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20C

Application Details

Western blotting:1:200-1000

Images



Gel: 6%SDS-PAGE

Lysate: 40 µg, Lane: K562 cell lysate

Primary antibody:CTR9?Antibody at dilution 1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 10MIN

Background

The protein encoded by this gene is a component of the PAF1 complex, which associates with RNA polymerase II and functions in transcriptional regulation and elongation. This complex also plays a role in the modification of histones.

Note: This product is for in vitro research use only