

CYB5R4 Antibody

Catalog No: #46999

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CYB5R4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total CYB5R4 protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human CYB5R4
Target Name	CYB5R4
Other Names	NCB5OR; cb5/cb5R; dJ676J13.1
Accession No.	Swiss-Prot#:Q7L1T6NCBI Gene ID:51167Gene Accssion:BC025380
Uniprot	Q7L1T6
GeneID	51167;
Concentration	0.5mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20C

Application Details

Immunofluorescence:1: 25-100

Images



The image is immunohistochemistry of paraffin-embedded Human liver cancer tissue using 46999(CYB5R4 Antibody) at dilution 1/20. (Original magnification: ?00)

Background

NADH-cytochrome b5 reductases are involved in desaturation and elongation of fatty acids, cholesterol biosynthesis, drug metabolism, and, in erythrocyte, methemoglobin reduction.

Note: This product is for in vitro research use only