

DEF8 Antibody

Catalog No: #47044

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DEF8 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total DEF8 protein.
Immunogen Type	protein
Immunogen Description	Full length fusion protein
Target Name	DEF8
Accession No.	Swiss-Prot#:Q6ZN54NCBI Gene ID:54849Gene Accssion:BC015482
Uniprot	Q6ZN54
GeneID	54849;
Concentration	0.7mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20C

Application Details

Immunofluorescence:1: 20-100

Images



The image is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using 47044(DEF8 Antibody) at dilution 1/25. (Original magnification: ?00)

Background

The function of this protein remains unknown.

Note: This product is for in vitro research use only