

## DENND5A Antibody

Catalog No: #47046

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	DENND5A Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total DENND5A protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human DENND5A
Target Name	DENND5A
Other Names	RAB6IP1
Accession No.	Swiss-Prot#:Q6IQ26NCBI Gene ID:23258Gene Accssion:BC009354
Uniprot	Q6IQ26
GeneID	23258;
Concentration	0.6mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20C

## Application Details

Immunofluorescence:1: 20-100

## Images



The image is immunohistochemistry of paraffin-embedded Human liver cancer tissue using 47046(DENND5A Antibody) at dilution 1/20. (Original magnification: ?00)

## Background

Guanine nucleotide exchange factor (GEF) which may activate RAB6A and RAB39A and/or RAB39B. Promotes the exchange of GDP to GTP, converting inactive GDP-bound Rab proteins into their active GTP-bound form.

---

Note: This product is for in vitro research use only