

DFNA5 Antibody

Catalog No: #47048

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DFNA5 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total DFNA5 protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human DFNA5
Target Name	DFNA5
Other Names	ICERE-1
Accession No.	Swiss-Prot#:O60443NCBI Gene ID:1687Gene Accssion:BC019689
Uniprot	O60443
GeneID	1687;
Concentration	0.6mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20C

Application Details

Immunofluorescence:1: 20-100

Images



The image is immunohistochemistry of paraffin-embedded Human tonsil tissue using 47048(DFNA5 Antibody) at dilution 1/20. (Original magnification: ?00)

Background

Hearing impairment is a heterogeneous condition with over 40 loci described. The protein encoded by this gene is expressed in fetal cochlea, however, its function is not known. Nonsyndromic hearing impairment is associated with a mutation in this gene. Three transcript variants encoding two different isoforms have been found for this gene.

Note: This product is for in vitro research use only