

## DHRS7 Antibody

Catalog No: #47052

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	DHRS7 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	WB, IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total DHRS7 protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human DHRS7
Target Name	DHRS7
Other Names	CGI-86; SDR34C1; retDSR4; retSDR4
Accession No.	Swiss-Prot#:Q9Y394NCBI Gene ID:51635Gene Accssion:BC000637
Uniprot	Q9Y394
GeneID	51635;
Calculated MW	38 kDa
Concentration	1mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20C


## Application Details

Western blotting:1:500-2000Immunofluorescence:1: 25-100


## Images



The image is immunohistochemistry of paraffin-embedded Human liver cancer tissue using 47052(DHRS7 Antibody) at dilution 1/30. (Original magnification: ?00)



The image is immunohistochemistry of paraffin-embedded Human prostate cancer tissue using 47052(DHRS7 Antibody) at dilution 1/30. (Original magnification: ?00)



Gel: 8%SDS-PAGE  
Lysate: 40 µg, Lane: 231 cell lysate  
Primary antibody:DHRS7 Antibody at dilution 1/400  
Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution  
Exposure time: 5 seconds

## Background

This gene encodes a member of the short-chain dehydrogenases/reductases (SDR) family, which has over 46,000 members. Members in this family are enzymes that metabolize many different compounds, such as steroid hormones, prostaglandins, retinoids, lipids and xenobiotics.?

Note: This product is for in vitro research use only