

DLL1 Antibody

Catalog No: #47062

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DLL1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	WB, IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total DLL1 protein.
Immunogen Type	peptide
Immunogen Description	Synthetic peptide of human DLL1
Target Name	DLL1
Other Names	DL1; Delta; DELTA1
Accession No.	Swiss-Prot#:O00548 NCBI Gene ID:28514Gene Accssion:NP_005609
Uniprot	O00548
GeneID	28514;
Calculated MW	78 kDa
Concentration	0.6mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20C


Application Details

Western blotting:1:200-1000Immunofluorescence:1: 20-100

Images



The image is immunohistochemistry of paraffin-embedded Human prostate cancer tissue using 47062(DLL1 Antibody) at dilution 1/20. (Original magnification: ?00)



Gel: 8%SDS-PAGE
Lysate: 40 µg, Lane: Rat heart tissue lysate
Primary antibody: DLL1 Antibody at dilution 1/250
Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution
Exposure time: 20 seconds

Background

DLL1 is a human homolog of the Notch Delta ligand and is a member of the delta/serrate/jagged family. It plays a role in mediating cell fate decisions during hematopoiesis. It may play a role in cell-to-cell communication.

Note: This product is for in vitro research use only