

ENO4 Antibody

Catalog No: #47109

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ENO4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	WB, IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total ENO4 protein.
Immunogen Type	peptide
Immunogen Description	Synthetic peptide of human ENO4
Target Name	ENO4
Other Names	C10orf134
Accession No.	Swiss-Prot#:A6NNW6NCBI Gene ID:387712Gene Accssion:NP_001229628
Uniprot	A6NNW6
GeneID	387712;
Calculated MW	69 kDa
Concentration	1.1mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20C


Application Details

Western blotting:1:500-2000Immunofluorescence:1: 20-100

Images



The image is immunohistochemistry of paraffin-embedded Human liver cancer tissue using 47109(ENO4 Antibody) at dilution 1/30. (Original magnification: ?00)



Gel: 8%SDS-PAGE
Lysate: 40 µg, Lane: Mouse lung tissue lysate
Primary antibody: ENO4 Antibody at dilution 1/450
Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution
Exposure time: 10 seconds

Background

May be required for sperm motility and function.

Note: This product is for in vitro research use only