

EXOSC3 Antibody

Catalog No: #47116

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	EXOSC3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total EXOSC3 protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human EXOSC3
Target Name	EXOSC3
Other Names	p10; PCH1B; RRP40; Rrp40p; CGI-102; hRrp-40; bA3J10.7
Accession No.	Swiss-Prot#:Q9NQT5NCBI Gene ID:51010Gene Accssion:BC008880
Uniprot	Q9NQT5
GeneID	51010;
Concentration	1mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20C


Application Details

Immunofluorescence:1: 25-100

Images



The image is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using 47116(EXOSC3 Antibody) at dilution 1/30. (Original magnification: ?00)



The image is immunohistochemistry of paraffin-embedded Human breast cancer tissue using 47116(EXOSC3 Antibody) at dilution 1/30. (Original magnification: ?00)

Background

This gene encodes a non-catalytic component of the human exosome, a complex with 3'-5' exoribonuclease activity that plays a role in numerous RNA processing and degradation activities. Related pseudogenes of this gene are found on chromosome 19 and 21. Alternatively spliced transcript variants encoding different isoforms have been described.

Note: This product is for in vitro research use only