

## FAM105A Antibody

Catalog No: #47118

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	FAM105A Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	WB, IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total FAM105A protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human FAM105A
Target Name	FAM105A
Other Names	NET20
Accession No.	Swiss-Prot#:Q9NUU6NCBI Gene ID:54491Gene Accsson:BC011524
Uniprot	Q9NUU6
GeneID	54491;
Calculated MW	42 kDa
Concentration	0.6mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20C


## Application Details

Western blotting:1:500-2000Immunofluorescence:1: 20-100

## Images



The image is immunohistochemistry of paraffin-embedded Human tonsil tissue using 47118(FAM105A Antibody) at dilution 1/20. (Original magnification: ?00)



Gel: 8%SDS-PAGE  
Lysate: 40 µg, Lane: Human placenta tissue lysate  
Primary antibody:FAM105A Antibody at dilution 1/400  
Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution  
Exposure time: 15 seconds

## Background

FAM105A belongs to the FAM105 family. The exact function of FAM105A remains unknown.

Note: This product is for in vitro research use only