

GPN2 Antibody

Catalog No: #47130

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GPN2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total GPN2 protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human GPN2
Target Name	GPN2
Other Names	ATPBD1B
Accession No.	Swiss-Prot#:Q9H9Y4NCBI Gene ID:54707Gene Accssion:BC008634
Uniprot	Q9H9Y4
GeneID	54707;
Concentration	0.7mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20C


Application Details

Immunofluorescence:1: 20-100

Images



The image is immunohistochemistry of paraffin-embedded Human cervical cancer tissue using 47130(GPN2 Antibody) at dilution 1/20. (Original magnification: ?00)



The image is immunohistochemistry of paraffin-embedded Human tonsil tissue using 47130(GPN2 Antibody) at dilution 1/20. (Original magnification: ?00)

Background

Small GTPase required for proper localization of RNA polymerase II and III (RNAPII and RNAPIII). May act at an RNAP assembly step prior to nuclear import.

Note: This product is for in vitro research use only