

GTF2F1 Antibody

Catalog No: #47131

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GTF2F1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total GTF2F1 protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human GTF2F1
Target Name	GTF2F1
Other Names	BTF4; RAP74; TF2F1; TFIIF
Accession No.	Swiss-Prot#:P35269NCBI Gene ID:2962Gene Accssion:BC013007
Uniprot	P35269
GeneID	2962;
Concentration	0.7mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20C

Application Details

Immunofluorescence:1: 20-100

Images



The image is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using 47131(GTF2F1 Antibody) at dilution 1/20. (Original magnification: ?00)

Background

TFIIF is a general transcription initiation factor that binds to RNA polymerase II and helps to recruit it to the initiation complex in collaboration with TFIIB. It promotes transcription elongation.

Note: This product is for in vitro research use only