

## MLLT3 Antibody

Catalog No: #47154

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	MLLT3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total MLLT3 protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human MLLT3
Target Name	MLLT3
Other Names	AF9; YEATS3
Accession No.	Swiss-Prot#:P42568NCBI Gene ID:4300Gene Accssion:BC036089
Uniprot	P42568
GeneID	4300;
Calculated MW	63 kDa
Concentration	0.1mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20C

## Application Details

Western blotting:1:200-1000

## Images



Gel: 8%SDS-PAGE  
 Lysate: 40 µg, Lane: 293T cell lysate  
 Primary antibody:MLLT3 Antibody at dilution 1/200  
 Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution  
 Exposure time: 5 minutes

## Background

Component of the super elongation complex (SEC), a complex required to increase the catalytic rate of RNA polymerase II transcription by suppressing transient pausing by the polymerase at multiple sites along the DNA.

---

Note: This product is for in vitro research use only