

NOLC1 Antibody

Catalog No: #47167

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NOLC1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total NOLC1 protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human NOLC1
Target Name	NOLC1
Other Names	P130; NOPP130; NOPP140; NS5ATP13
Accession No.	Swiss-Prot#:Q14978NCBI Gene ID:9221Gene Accssion:BC006769
Uniprot	Q14978
GeneID	9221;
Concentration	0.9mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20C

Application Details

Immunofluorescence:1: 25-100

Images



The image is immunohistochemistry of paraffin-embedded Human tonsil tissue using 47167(NOLC1 Antibody) at dilution 1/25. (Original magnification: ?00)

Background

Related to nucleogenesis, may play a role in the maintenance of the fundamental structure of the fibrillar center and dense fibrillar component in the nucleolus. It has intrinsic GTPase and ATPase activities. May play an important role in transcription catalyzed by RNA polymerase I.

Note: This product is for in vitro research use only