

TRMT1 Antibody

Catalog No: #47226

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TRMT1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total TRMT1 protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human TRMT1
Target Name	TRMT1
Other Names	TRM1
Accession No.	Swiss-Prot#:Q9NXH9NCBI Gene ID:55621Gene Accssion:BC002492
Uniprot	Q9NXH9
GeneID	55621;
Concentration	0.8mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20C


Application Details

Immunofluorescence:1: 20-100

Images



The image is immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using 47226(TRMT1 Antibody) at dilution 1/25. (Original magnification: ?00)



The image is immunohistochemistry of paraffin-embedded Human liver cancer tissue using 47226(TRMT1 Antibody) at dilution 1/25. (Original magnification: 200)

Background

Dimethylates a single guanine residue at position 26 of most tRNAs using S-adenosyl-L-methionine as donor of the methyl groups.

Note: This product is for in vitro research use only