

TSSK3 Antibody

Catalog No: #47228

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TSSK3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total TSSK3 protein.
Immunogen Type	protein
Immunogen Description	Full length fusion protein
Target Name	TSSK3
Other Names	TSK3; SPOGA3; STK22C; STK22D
Accession No.	Swiss-Prot#:Q96PN8NCBI Gene ID:81629Gene Accssion:BC035354
Uniprot	Q96PN8
GeneID	81629;
Concentration	0.9mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20C


Application Details

Immunofluorescence:1: 25-100

Images



The image is immunohistochemistry of paraffin-embedded Human liver cancer tissue using 47228(TSSK3 Antibody) at dilution 1/25. (Original magnification: ?00)



The image is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using 47228(TSSK3 Antibody) at dilution 1/25. (Original magnification: ?00)

Background

This gene encodes a kinase expressed exclusively in the testis that is thought to play a role in either germ cell differentiation or mature sperm function.

Note: This product is for in vitro research use only