

TUSC2 Antibody

Catalog No: #40349

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TUSC2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total TUSC2 protein.
Immunogen Type	Protein
Immunogen Description	Full length fusion protein
Target Name	TUSC2
Other Names	PAP; FUS1; PDAP2; C3orf11
Accession No.	Swiss-Prot:O75896Gene Accssion:BC023976
Uniprot	O75896
GeneID	11334;
Concentration	2.9mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry:1:30-1:150

Images



Immunohistochemical analysis of paraffin-embedded Human liver cancer tissue using #40349 at dilution 1/50.

Background

This gene is a highly conserved lung cancer candidate gene. No other information about this gene is currently available.

Note: This product is for in vitro research use only