

UACA Antibody

Catalog No: #40351

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	UACA Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total UACA protein.
Immunogen Type	Protein
Immunogen Description	Fusion protein of human uveal autoantigen with coiled-coil domains and ankyrin repeats
Target Name	UACA
Other Names	NUCLING
Accession No.	Swiss-Prot:Q9BZF9Gene Accssion:BC113407
Uniprot	Q9BZF9
GeneID	55075;
Concentration	2.7mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry:1:30-1:150

Images



Immunohistochemical analysis of paraffin-embedded Human liver cancer tissue using #40351 at dilution 1/45.

Background

Regulates APAF1 expression and plays an important role in the regulation of stress-induced apoptosis. Promotes apoptosis by regulating three pathways, apoptosome up-regulation, LGALS3/galectin-3 down-regulation and NF-kappa-B inactivation. Regulates the redistribution of APAF1 into the nucleus after proapoptotic stress. Down-regulates the expression of LGALS3 by inhibiting NFKB1 (By similarity).

Note: This product is for in vitro research use only