

FOXG1 Antibody

Catalog No: #40353

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	FOXG1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous levels of total FOXG1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide corresponding to residues near the C terminal of human forkhead box G1
Target Name	FOXG1
Other Names	BF1; BF2; QIN; FKH2; HBF2; HFK1; HFK2; HFK3; KHL2; FHKL3; FKHL1; FKHL2; FKHL3; FKHL4; HBF-1; HBF-2; HBF-3; FOXG1A; FOXG1B; FOXG1C; HBF-G2
Accession No.	Swiss-Prot:P55316Gene Accssion:NP_005240
Uniprot	P55316
GeneID	2290;
SDS-PAGE MW	52KD
Concentration	2.1mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:200-1:1000

Images



Gel: 8%SDS-PAGE
 Lysates (from left to right): Jurkat, K562 and 293T cell
 Amount of lysate: 40ug per lane
 Primary antibody: 1/200 dilution
 Secondary antibody dilution: 1/8000
 Exposure time: 5 minutes

Background

This locus encodes a member of the forked-head transcription factor family. The encoded protein, which functions as a repressor, may play a role in

brain development. Mutations at this locus have been associated with Rett syndrome.?

Note: This product is for in vitro research use only