

TNFSF4 Antibody

Catalog No: #40355

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TNFSF4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total TNFSF4 protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide of human tumor necrosis factor (ligand) superfamily, member 4
Target Name	TNFSF4
Other Names	GP34; CD252; OX40L; TXGP1; CD134L; OX-40L
Accession No.	Swiss-Prot:P23510Gene Accssion:NP_003317
Uniprot	P23510
GeneID	7292;
Concentration	1.4mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry:1:10-1:50

Images



Immunohistochemical analysis of paraffin-embedded Human liver cancer tissue using #40355 at dilution 1/20.



Immunohistochemical analysis of paraffin-embedded Human brain tissue using #40355 at dilution 1/20.

Background

This gene encodes a cytokine of the tumor necrosis factor (TNF) ligand family. The encoded protein functions in T cell antigen-presenting cell (APC) interactions and mediates adhesion of activated T cells to endothelial cells. Polymorphisms in this gene have been associated with Sjogren's syndrome and systemic lupus erythematosus. Alternative splicing results in multiple transcript variants.

Note: This product is for in vitro research use only