

AEBP1 Antibody

Catalog No: #40358

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	AEBP1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total AEBP1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide of human AE binding protein 1
Target Name	AEBP1
Other Names	ACLP
Accession No.	Swiss-Prot:Q8IUX7Gene Accssion:NP_001120
Uniprot	Q8IUX7
GeneID	165;
Concentration	1.5mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry:1:25-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human esophagus cancer tissue using #40358 at dilution 1/25.



Immunohistochemical analysis of paraffin-embedded Human thyroid cancer tissue using #40358 at dilution 1/25.

Background

This gene encodes a member of carboxypeptidase A protein family. The encoded protein may function as a transcriptional repressor and play a role in adipogenesis and smooth muscle cell differentiation. Studies in mice suggest that this gene functions in wound healing and abdominal wall development. Overexpression of this gene is associated with glioblastoma.

Note: This product is for in vitro research use only