

TPST2 Antibody

Catalog No: #40389

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TPST2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total TPST2 protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide corresponding to residues near the N terminal of human tyrosylprotein sulfotransferase 2
Target Name	TPST2
Other Names	TANGO13B
Accession No.	Swiss-Prot:O60704Gene Accssion:NP_003586
Uniprot	O60704
GeneID	8459;
Concentration	1.3mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Colorectal cancer tissue using #40389 at dilution 1/100.



Immunohistochemical analysis of paraffin-embedded Human Lung cancer tissue using #40389 at dilution 1/100.

Background

The protein encoded by this gene catalyzes the O-sulfation of tyrosine residues within acidic regions of proteins. The encoded protein is a type II integral membrane protein found in the Golgi body. Two transcript variants encoding the same protein have been found for this gene.

Note: This product is for in vitro research use only