

Mouse IGF-2 ELISA Kit

Catalog No: #EK5155



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse IGF-2 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Ms
Target Name	Igf2
Other Names	Insulin-like growth factor II; IGF-II; Multiplication-stimulating polypeptide; Insulin-like growth factor II; Preptin; Igf2; Igf-2;
Accession No.	P09535
Uniprot	P09535
GeneID	16002;

Application Details

Detect Range: 62.5 pg/mL - 4000 pg/mL

Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul

Assay Time: 3-5h

Detection wavelength: 450 nm

Note: This product is for in vitro research use only