

PEPC (Phospho-Ser15) Polyclonal Antibody

Catalog No: #12178

Package Size: #12178-1 50ul #12178-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PEPC (Phospho-Ser15) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from antiserum by affinity-chromatography
Applications	WB
Species Reactivity	Plant
Specificity	The antibody detects endogenous PEPC (Phospho Ser15) protein.
Immunogen Type	peptide
Immunogen Description	Phospho peptide at Ser15 of PEPC
Target Name	PEPC
Modification	Phospho
Other Names	PGC; Gastricsin; Pepsinogen C
Accession No.	Swiss-Prot: P20142NCBI Gene ID: 5225
Uniprot	P20142
GeneID	5225;
SDS-PAGE MW	105-110kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:1000-2000

Images



Western blot analysis of Maize 1) Dark, 2) Light, diluted at 1:2000.

Note: This product is for in vitro research use only